



**Sensitivity  
Flexibility  
Experience**

# A New Generation of Evaporative Light-Scattering Detectors for Liquid Chromatography: Universality, High Performance and Robustness in Pharmaceutical Analysis - An Application Review in HPLC and U-HPLC

## Abstract


Among the detectors available in Liquid Chromatography (LC), Evaporative Light-Scattering Detector (ELSD) became in recent years a well established instrument thanks to several theoretical studies based on fundamental investigations and numerous applications provided during the last thirty years. Indeed, ELSD is considered as a nearly Universal, powerful, reliable and cost-effective technique, and is ideally appropriate in Pharmaceutical industry for a great variety of LC applications containing chromophoric and non-chromophoric compounds.

Today, an ELSD model based on a recent and unique concept is proposed which offers a genuine and efficient Low-Temperature technology (LT-ELSD™) combined to an innovative detection chamber. The overall design of this ultimate detector results in a significant increase of sensitivity providing typical limits of detection down to the very low nanogram levels for non-volatile and semi-volatile compounds. It provides an improved overall direct linearity with correlation coefficients over 0.99, consistent responses independent of the analytes chemical structure and an extended dynamic range exceeding the four orders of magnitude (from low ng to high µg levels on column). Also, this model is optimized for the recent U-HPLC technique giving peak widths of less than 1 second.

To show the strength and the versatility of this ELSD model, several relevant LC applications in Pharmaceutical analysis are developed in this work. These applications use the most recent LC media, such as multi-mode, HILIC and sub-two-micron or fused-core particle phases, allowing outstanding separations and simultaneous analyses of a wide range of chemical and biochemical compounds.

The topics presented here are focused on:


- Sensitivity and robustness in LC on the full pH range.
- Rapid and simultaneous separation of API, counterions and impurities.
- Response consistency compared to UV detectors at several wavelengths.
- Relevance in the determination of both chromophoric and non-chromophoric solutes in the analysis of natural products or TCM, such as Ginkgo Biloba.
- Simplified and cost-effective alternative in the analysis of non-chromophoric compounds such as aminoglycoside antibiotics, thus avoiding any tedious derivatization step, ion pair reagents and specific detectors.



**HPLC / LT-ELSD System**

## I - Sensitivity

Application: Generic HPLC/LT-ELSD Analysis of Several Lipid Compounds



	RT	LOD (S/N=3)
	Minutes	ng (o.c.)
1 - Lauric acid	4.87	16.2*
2 - Linolenic acid	7.17	4.1
3 - Myristic acid	7.58	1.6
4 - Retinol (Vit A)	8.10	3.6
5 - Linoleic acid	9.43	5.1
6 - Monolein	10.21	4.8
7 - Palmitic acid	11.43	0.8
8 - Oleic acid	12.35	5.7
9 - Hexadecanol	12.88	2.1
10 - Stearic acid	15.77	0.5
11 - Octadecanol	17.32	0.5
12 - Eicosanol	21.63	0.7
13 - Cholesterol	23.80	1.3
14 - Docosanol	25.57	0.9
15 - α-Tocopherol (Vit E)	25.80	3.8
16 - Vitamin K	29.20	3.8
17 - Squatene	32.54	2.4
18 - Diolefin	34.13	2.3
19 - Trilaurin	36.50	2.1
20 - Trilinolenin	38.90	2.5
21 - Trimyristin	40.97	1.7
22 - Coenzyme Q10	41.09	1.8
23 - Trilinolein	41.73	1.9
24 - Tripalmitin	44.09	1.7
25 - Triolein	44.29	1.1

\* Semi-volatile compound

Chromatogram of the Simultaneous HPLC/LT-ELSD Analysis of Fatty Acids, Fatty Alcohols, Fat-Soluble Vitamins, Mono-, Di- and TriGlycerides and Related Compounds.


Standard mixture: 25 Compounds (see Table beside)  
Injection volume: 2µL  
Column: Hypersil GOLD (1.9µm, 2.1 x 200mm), 60°C  
Flowrate: 0.3mL/min  
Eluent: A: MeOH/ACN/CH<sub>2</sub>Cl<sub>2</sub>/formic acid (500:300:198:2) - B: MeOH/acetone/formic acid (998:400:2)  
Gradient: 0-3 minutes: 100%A, 3-43 minutes: from 100%A to 100%B

Detector: SEDEX 90LT, 28°C, 3.5Bar

The results show very high sensitivities obtained in a real HPLC/LT-ELSD application. Obtained LODs are much below 10ng on column for all compounds (except for Lauric acid which is characterized by a high vapor pressure), and even at the Picogram Levels for some other semi-volatile compounds belonging to the groups of fatty alcohols and fatty acids.

## II - API and their Counterions

Application: Global HPLC/LT-ELSD Method for the Simultaneous Analysis of Polar and Non-Polar, Neutral, Acidic and Basic Pharmaceutical Drugs and their Respective Counterions



Multimodal Stationary Phase HPLC/LT-ELSD Chromatogram of the Simultaneous Analysis of Polar and Non-Polar, Neutral, Acidic and Basic Pharmaceutical Drugs and their Counterions.

Standard mixture: 1- Acetaminophen, 2- Sodium, 3- Potassium, 4- Hydrocortisone, 5- Procainamide, 6- Chloride, 7- Nitrate, 8- Micronazole, 9- Lorazepam, 10- Dichlorofenac (500ppm each API)  
Injection volume: 2µL  
Column: Acclaim Trinity P1 (3µm, 2.1 x 150mm), 30°C  
Flowrate: 0.35mL/min  
Eluent: 80% Ammonium acetate 20mM, pH5 + 20% ACN (A) / 30% Ammonium formate, 20mM, pH3 + 70% ACN (B)  
Gradient: 0-2 minutes: 95%B, 2-17 minutes: from 95%B to 100%B


Detector: SEDEX 90LT, 40°C, 3.5Bar

This chromatogram shows that the combination of an efficient multimodal stationary phase and a single Universal ELSD detector allows the quick and easy simultaneous determination of a wide range of compounds characterized by different polarities, including inorganic anions and cations.

**Eric VERETTE, Ph.D.**  
SEDERE S.A.S., France

## III - Impurity Assessment

Application: Simultaneous HPLC/LT-ELSD Analysis of Imipramine, its Counterion and an Impurity



Chromatograms of the Simultaneous HPLC/LT-ELSD Analysis of Imipramine and its Counterion, with and without an Impurity (Sodium, 5ppm)


Standard mixture: 1- Imipramine (API: 10 000ppm), 2- Cl<sup>-</sup> (Counterion)  
Injection volume: 2µL (20µg API)  
Column: Acclaim Trinity P1 (3µm, 2.1 x 150mm), 35°C  
Flowrate: 0.5mL/min  
Eluent: Ammonium acetate 50mM, pH5 / ACN (60:40)

Detector: SEDEX 90LT, 40°C, at 0 minute: Gain 10, at 2 minutes: Gain 4, 3.5Bar

In this example SEDEX 90LT detects an impurity (sodium) at a level of 0.05% of the major compound (Imipramine). With this detector, there are 12 gains available, which means that this percentage could typically go down to less than 0.01% (more than 4 orders of magnitude for the overall dynamic range).

## IV - Response Consistency

Application: Fast HPLC/LT-ELSD/DAD Analysis of Non-Volatile Compounds with Different Chemical Structures




Chromatograms of the HPLC/LT-ELSD/DAD Analysis of 5-Fluorocytosine, Theophylline and Acetaminophen

Standard Mixture: 1- 5-Fluorocytosine, 2- Theophylline and 3- Acetaminophen (500ppm each)  
Injection volume: 2µL  
Column: Halo C18 (2.7µm, 2.1 x 150mm), 30°C  
Flowrate: 0.5mL/min  
Eluent: H<sub>2</sub>O / ACN (85:15)

Detector: SEDEX 90LT, 50°C, 3.5Bar

## V - Robustness at High pH

Application: Simultaneous U-HPLC/LT-ELSD Analysis of Several Beta Blockers and Tricyclic Antidepressants at pH:11



Chromatogram of the Simultaneous U-HPLC/LT-ELSD Analysis of Several Beta Blockers and Tricyclic Antidepressants


Standard Mixture: 1- Atenolol, 2- Pindolol, 3- Acebutolol, 4- Metoprolol, 5- Propranolol, 6- Nortriptyline, 7- Imipramine, 8- Amitriptyline (500ppm each)  
Injection volume: 1µL  
Column: Zorbax Extend C18 (1.8µm, 2.1 x 50mm), 40°C  
Flowrate: 0.3mL/min  
Eluent: H<sub>2</sub>O + Triethylamine 20mM, pH:11 (A) / Methanol (B)  
Gradient: 0-0.5 minute: 35%B, 0.5-4 minutes: from 35%B to 95%B, 4-8 minutes: 95%B

Detector: SEDEX 90LT, 40°C, 3.5Bar

This example shows that SEDEX 90LT is not affected by very basic buffers such as 20mM triethylamine (pH:11), and provides a nice and flat baseline with no drift all along the gradient.

## VI - Alternative Solution to Derivatization and Ion-Pair Reagents

Application: Direct HILIC/LT-ELSD Analysis of Aminoglycoside Antibiotics



Chromatogram of 6 Aminoglycoside Antibiotics by HILIC/LT-ELSD


Standard mixture: 1- Streptomycin, 2- Amikacin, 3- Kanamycin, 4- Paromomycin, 5- Tobramycin, 6- Neomycin (1000ppm each)  
Injection volume: 2µL  
Column: Ascentis Express HILIC (2.7µm, 2.1 x 150mm), 30°C  
Flowrate: 0.6mL/min  
Eluent: Ammonium acetate 150mM, pH5 (A) / ACN (B)  
Gradient: 0-1 minute: 60%B, 1-5 minutes: from 60%B to 5%B, 5-7 minutes: 5%B

Detector: SEDEX 90LT, 60°C, 3.5Bar

SEDEX 90LT allows the direct determination of aminoglycoside antibiotics. With these chromatography conditions, Gentamicin congeners (not shown here) elute after 7 minutes.

## VII - Analysis of Natural Products

Application: Simultaneous U-HPLC/LT-ELSD Analysis of Terpenic Lactones and Flavonoids Contained in Ginkgo Biloba



Chromatogram of 4 Terpenic Lactones and 3 Flavonoids by U-HPLC/LT-ELSD

Standard mixture: 1- Bilobalide, 2- Ginkgolide C, 3- Ginkgolide A, 4- Ginkgolide B, 5- Quercetin, 6- Isohannitinin, 7- Kaempferol (250ppm each)  
Injection volume: 1µL  
Column: Hypersil GOLD (1.9µm, 2.1 x 50mm), 30°C  
Flowrate: 0.6mL/min  
Eluent: H<sub>2</sub>O + 0.1% formic acid / Acetone + 0.1% formic acid  
Gradient: 0-0.5 minute: 5%B, 0.5-4 minutes: from 5%B to 50%B, 4-6 minutes: 50%B

Detector: SEDEX 90LT, 50°C, 3.5Bar

SEDEX 90LT allows the simultaneous determination of both chromophoric (flavonoids) and non-chromophoric (terpenic lactones) compounds, which shows the great advantage of this single, easy to operate and Universal detection mode. Elsewhere, it allows the use of acetone which is less toxic and cheaper than many other organic eluents (such as acetonitrile), and which possesses excellent physical and chromatography properties.

## Conclusion

The applications developed here clearly show the advantages of the new SEDEX 90LT ELSD and particularly in regards to:

- Sensitivity with low nanogram and even sub-nanogram levels.
- Simultaneous analysis of both chromophoric and non-chromophoric solutes, using just a single Universal detector.
- Wide dynamic range allowing a sensitive impurity assessment.
- Very small response variation between compounds, compared to UV detectors.
- Use of acetone in the mobile phase, which cannot be selected with UV detectors due to the high cutoff.

This work also demonstrates the significant advancement of the new SEDEX Evaporative Light-Scattering Detector resulting from the combination of an efficient and genuine Low Temperature technology and an innovative detection device based on a high-performance laser. These outstanding new features offer now to the analyst a Universal, powerful, versatile and cost-effective solution to their separation and quantification challenges in Pharmaceutical area.