


## Rapid production and size assessment of Embryoid Bodies ( EBs )



4 Plate type




Single plate type



### Simplify formation of EBs with GravityPLUS™ hanging drop plates and Cell<sup>3</sup>iMager analysis

- Produce consistently sized EBs in a 96-well automation-compatible format
- Eliminate suspension culture and tedious manual selection of EBs with simple transfer to culture plates
- Monitor EB growth and differentiation with rapid, label-free imaging



Size analysis of embryoid body at day 5 of culture using the Cell<sup>3</sup>iMager

### A no hassle solution for production of consistently sized EBs

The establishment of 3-dimensional embryoid bodies (EBs) from pluripotent stem cells holds promise as tool for regenerative medicine, in vitro toxicity testing, and as a model for studying embryogenesis.

Large-scale production of EBs for use in such purposes can be accomplished by seeding single ESCs in rotating suspension culture, a process resulting in large numbers of EBs, but limited in its ability to reliably produce EBs of consistent size. Alternatively, more consistently sized EBs can be formed using the hanging drop method, in which single-cell suspensions are seeded and grown on petri dish lids, transferred to suspension culture, then hand-picked for analysis in high-throughput culture plates.

However this tedious, labor-intensive process can also result in size variations resulting from fusion of EBs during the suspension culture phase.

InSphero's patented GravityPLUS™ hanging drop platform provides a simple, automation-compatible solution to the production of consistently sized assay-ready EBs.

The SureDrop™ inlet allows hassle-free formation of EBs in hanging drop by simple top-loading of cells in suspension.

The stable 96-well format increases throughput, reduce waste, and eliminates the need for suspension culture and cumbersome hand-picking of similarly sized EBs. Subsequent analysis with SCREEN Cell<sup>3</sup>iMager provides label-free monitoring of EB size and morphology over days in culture, ensuring consistency of EB production and as an endpoint for growth and differentiation.

#### Applications




- Embryotoxicity / teratogenicity screening
- Differentiation pathway analysis
- Quality control

# Rapid production of Embryoid Bodies Never flip your lid again!




## Simplify formation of EBs with GravityPULS™ hanging drop plates

ESD3[D3](ATCC® CRL1934™) were expanded on gelatin coated flasks in serum containing stem cell medium supplemented with LIF. Subsequently, ESD3 cells were trypsinized and harvested to re-seed in the 96-well GravityPLUS™ hanging drop platform at densities of 250, 500, 750 and 1000 cells/well (24 replicates per density). Cells were allowed to re-aggregate to form EBs in a serum supplemented cell differentiation medium.



## Assessment of EB size in GravityTRAP™ plates using the Cell<sup>3</sup>iMager

After 3 days, developed EBs were transferred into the GravityTRAP™ non-adherent flat bottom plate by addition on fresh culture medium through the SureDrop™ inlet. EBs were monitored for growth at 3 time points (days 3, 4 and 5) with an optical read-out using the Cell<sup>3</sup>iMager (SCREEN Holdings Co., Ltd. Japan) at a resolution of 4800 dpi. Captured images were used to calculate average EB area (top), diameter (bottom) and optical density (data not shown) applying an algorithm specific for EB size assessment.



Top: overview of EB formation in the GravityPLUS™ hanging drop platform and subsequent transfer to GravityTRAP™ plates for size analysis on Cell<sup>3</sup>iMager. Bottom: representative images depicting EBs following transfer and culture in 96 well-GravityTRAP™ plate (day 3 = day of transfer, well diameter = 1mm)

**InSphero AG**  
Schlieren, Switzerland +41-44-515049-0

**InSphero USA Inc.**  
Cambridge, MA, USA +1-800-779-7558

**InSphero Europe GmbH**  
Waldshut, Germany +49-7751-304966-5

**Authorized Distributors**  
[www.insphero.com/distributors](http://www.insphero.com/distributors)  
[info@insphero.com](mailto:info@insphero.com) [www.insphero.com](http://www.insphero.com)

InSphero is ISO 9001:2008 certified  
All rights reserved, 2013 InSphero AG. 3D InSight™ and GravityPLUS™ are trademarks of InSphero AG. For life science research only. Not for use in diagnostic procedures.

**SCREEN Holdings Co., Ltd.**  
Distributed by SciPro  
[www.scipro.com](http://www.scipro.com)  
[info@scipro.com](mailto:info@scipro.com)  
800-268-1956  
[www.screen.co.jp](http://www.screen.co.jp)